Quantitation of Protein by Colormetric Assays

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- 1) Dissolve a frozen aliquot of 1 mg/ml BGG for a protein standard (Prepare 5 dilutions of a protein standard containing 0.2 mg/ml to 1.5 mg/ml in the same buffer as the sample for the best results).
- 2) Prepare reagent A' by mixing 300 µl of reagent A and 6 µl of reagent S (or 25 µl x sample and standard No.).
- 3) Pipet ≤5 µl of standards and samples into microtiter plate wells.
- 4) Add 22.5 µl of reagent A' into each well.
- 5) Add 180 µl reagent B into each well, and gently agitate the plate to mix the reagents.
- 6) After 15 min, measure absorbance at 655 nm with a Microplate Reader.

Table 1. Reagents Compatible with the Bio-Rad Protein Assay When Using the Standard Procedure.*

Acetate, 0.6 M	Acetone	Adenosine, 1 mM	AminoAcids
Ammonium sulfate, 1.0 M	Ampholytes, 0.5%	Acid pH	ATP, 1 mM
Barbital	BES, 2.5 M	Boric acid	Cacodylate-Tris, 0.1 M
CDTA, 0.05 M	Citrate, 0.05 M	Deoxycholate, 0.1%	Dithiothreitol, 1 M
DNA, 1 mg/ml	EDTA, 0.1 M	EGTA, 0.05 M	Ethanol

Eagle's MEM	Earle's salt solution	Formic acid, 1.0 M	Fructose
Glucose	Glutathione	Glycerol, 99%	Glycine, 0.1 M
Guanidine-HCl	Hank's salt solution	HEPES buffer, 0.1 M	KCl, 1.0 M
Malic acid, 0.2 M	MgCl ₂ , 1.0 M	Mercaptoethanol, 1.0 M	
MES, 0.7 M	Methanol	MOPS, 0.2 M	NaCl, 5 M
NAD, 1 mM	NaSCN, 3 M	Peptones	Phenol, 5%
Phosphate, 1.0 M	PIPES, 0.5 M	Polyadenylic acid, 1 mM	
Polypeptides (MW<3000)	Pyrophosphate, 0.2 M	rRNA, 0.25 mg/ml	tRNA, 0.4 mg/ml
total RNA, 0.30 mg/ml	SDS, 0.1%	Sodium phosphate	
Streptomycin sulfate, 20%	Triton X-100, 0.1%	Tricine	Tyrosine, 1 mM
Thymidine, 1 mM	Tris, 2.0 M	Urea, 6 M	Vitamins

* Interference may be caused by chemical-protein and/or chemical-dye interactions. Table 1 lists those chemical reagents not directly affecting the development of dye color. Since every protein-chemical reagent combination has not been assayed, it is possible that some of the listed reagents produce interference in combination with certain proteins. However, with respect to proteins such as bovine albumin and globulin, the above listed reagents show little or no interference.

Table 2. Reagents compatible with the Bio-Rad DC Protein Assay.

In some cases, the presence of one or more of these substances will effect a change in the response of the protein to the assay reagents; therefore, the standard should always be prepared in the same buffer as the sample.

10% SDS	1% CHAPS	2% NP-40	1% Triton X-100
1% CHAPSO	1% Thesit	1% Tween 20	1% Octyl glucoside
1% Brij-35	0.2% C12E8 *	0.1 M Tris, pH 8	0.5 M NaOH
0.5 M HCl	0.5 M (NH ₄) ₂ SO ₄	0.025 M EDTA	0.05 M CaCl ₂
0.4 M Guanidine HCl	4 M Urea	0.05% Sodium azide	1 mM DTT (dithiothreitol)

Note: The DC Protein Assay is incompatible with 2-mercaptoethanol (BME).

*octaethyleneglycol dodecyl ether