Synchronization/ 細胞同調

Rumi Miyamoto 020728 (+ Masuda)

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Preculture

Transfer a single colony from a streak plate to 3 ml YEA.
Incubate overnight at 30°C with shaking.

Culture

Determine cell number by counting.
Inoculate 60 ml YEA.¹⁾
Grow at 30°C with shaking for ~16 hr.

G1 arrest in EMM2-N culture

Determine cell number by counting.²⁾
Wash twice with EMM2-N.
Adjust cell concentration to 2.0×10^6 cells/ml by diluting with EMM2-N.
Incubate overnight (about 4 generation time) at 25°C with shaking.³⁾

Release

Resuspend cells in YEA at 2.0×10^6 cells/ml, and incubate at 33°C.

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Synchronization by HU treatment

Determine cell number by counting.²⁾

Adjust cell concentration to 2.0×10^6 cells/ml by diluting with YEA.

Incubate with 15 mM HU³⁾ in YEA for 3 hr.

Determine cell number by counting.⁴⁾

Collect cells and wash twice with YEA.

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Preculture

Transfer a single colony from a streak plate to 3 ml YEA.

Incubate overnight at 30 °C with shaking.

Culture

Innoculate 250 ml YEA \times 2 flasks.¹⁾

Grow at 30°C with shaking for ~16 hr.

Selection Synchrony

(Make 10-40% lactose gradient in 40 ml YEA.)²⁾

Determine cell number by counting.

(use cells at mid-exponential phase: $< 1 \times 10^7$ cells/ml)

Centrifuge at 3,500 rpm for 5 min.

Resuspend in YEA.³⁾ (total cell suspension: ≤ 3 ml)

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Apply cells³⁾ on lactose gradient medium.

Centrifuge at 2,000 rpm for 2 min with a swing rotor.³⁾

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Resuspend cells in YEA at 2×10^6 cells/ml, and incubate at 33°C.
(early/mid-log phase の細胞濃度でリリース)

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Preculture

Inoculate 10 ml YEA with a single colony from a YEA plate.

Incubate overnight at 30°C with shaking.

Culture

Determine cell number by counting.¹⁾

Inoculate 700 ml of YEA.²⁾

Incubate overnight at 30°C with shaking.

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Determine cell number by counting.⁴⁾

Collect cells by centrifuging at 4,000 rpm for 5 min.

Resuspend pellet in ~60 ml of YEA at 3.0×10^8 cells/ml.

Load 30 ml cell suspension ($3.0 \times 10^8 \times 30$ cells) into a 40 ml sample chamber.

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D) Synchronization by temperature block

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Cell : HR348-1 *h⁻ leu1 cdc25-22*

Preculture

Transfer a single colony from a streak plate to 3 ml YEA.

Incubate overnight at 25°C with shaking.

Culture

Determine cell number .

Innoculate 100 ml YEA.¹⁾

Grow at 26°C with shaking for ~16 hr.

Synchronization by temperature block

Determine cell number.

Adjust cell concentration to 3.0×10^6 cells/ml by diluting with YEA.

Incubate at 36°C for 4 hr.²⁾

Cool to 25°C by incubating the flask on ice with shaking.³⁾

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